



Installation Services

Backed by over a century of technical expertise, STERIS is your global leader and trusted partner for installing your equipment.

We'll work with you to coordinate every aspect of your equipment installation, including transportation, delivery and any facility-specific requirements. With both full installation packages and stand-alone services, STERIS has the expertise and resources to make your installation process quick and easy.

Full Installation Service

Our installation services range from stand-alone services to a full installation package. Our full package includes the necessary labor and materials to move your equipment, connect to the required utilities and verify equipment functionality. Individual services are also available and can be customized or combined to meet your facility's needs.

SET IN PLACE

Includes necessary labor and materials to:

- Uncrate, inspect and assemble equipment
- Relocate to the specified place your facility

FINAL CONNECTIONS

Includes labor and materials to make utility connections:

- Air
- Steam
- Water
- Drain
- Electric

START UP

Includes testing and inspection to verify:

- Equipment is properly installed
- Installation and utility connections meet specifications
- Operational cycles run correctly



Other Installation Services

STERIS offers a variety of additional services available to help meet regional and other specialized installation needs. We provide deinstallation, installation of seismic hardware and bioseals, coverage of local permits or inspections, and supervision of third-party installers to ensure your investment is protected and installed properly.

Contact a STERIS representative for help choosing the right installation services for your facility.