



## Heat Distribution Testing Service

Ensure your sterilizer delivers consistent, accurate results. This service is ideal for facilities with STERIS sterilization equipment that require validation of new equipment or verification of existing systems.

This testing identifies temperature deviations and slow-to-heat areas within the empty chamber to optimize cycle efficiency, validate performance and ensure reliable operation.

### Compatible STERIS Equipment Offerings:

- AMSCO™ Small & Medium Sterilizers
- AMSCO™ Bulk Sterilizers
- Finn-Aqua™ Sterilizers

### Your Service Includes



**APPROVED  
PROTOCOL**



**TRACEABLE  
CALIBRATION**



**PROTOCOL  
EXECUTION**



**FINAL  
REPORT**

All company and product names are trademarks of STERIS, its affiliates or related companies, unless otherwise noted. ©2025 STERIS