

Why Use An Acid Detergent?



Maintain Surface Finish

A rough surface finish is susceptible to potential corrosion, bacterial attachment and biofilm formation on equipment components and associated piping.



Improve Cleanability

Rouge prevents cleaning agent contact with underlying residues or soils, making surfaces more difficult to clean.



Prevent Microbial Excursions

A poor cleaning regime can lead to multi-species Biofilm:

Ralstonia pickettii

Pseudomonas fluorescens

Bacille non revivifiable

Burkholderia cepacia



Avoid Product Contamination

The presence of rouge or bacterial contamination due to rouge can impact product production or possibly adulterate the product.

All company and product names are trademarks of STERIS, its affiliates or related companies, unless otherwise noted. ©2025 STERIS